

## Background

Both Glucagon-Like Peptide-1 Receptor Agonists (GLP-1)<sup>1</sup> and continuous glucose monitors (CGM)<sup>2</sup> provide glycemic benefits for type 2 diabetes (T2D). However, it is unknown if adding CGM to GLP-1 therapy further improves HbA1c.

## Objective

In a real-world analysis, evaluate changes in HbA1c after acquisition of FreeStyle Libre (FSL) in adults with T2D on prior GLP-1 therapy.

## Methods

**Data Source:** Optum's de-identified Market Clarity Data, a linked electronic health records-claims database (2017-2022) based in the United States.

### Inclusion Criteria:


- FSL CGM acquisition
- T2D diagnosis and baseline HbA1c  $\geq 8\%$
- 18+ years of age
- 180 days of continuous Rx and medical coverage prior to FSL acquisition
- At least one GLP-1 acquisition within 180 days before FSL

### Exclusion Criteria:

- History of non-FSL CGM
- Gestational or type 1 diabetes

<sup>1</sup>Sibony RW, et al. Drug Therapies for Diabetes. Int J Mol Sci. 2023;24(24):17147. <sup>2</sup>Seidu S, et al. Efficacy and Safety of Continuous Glucose Monitoring and Intermittently Scanned Continuous Glucose Monitoring in Patients With Type 2 Diabetes: A Systematic Review and Meta-analysis of Interventional Evidence. Diabetes Care. 2024;47(1):169-179.

## HbA1c Measurements:



### Baseline HbA1c:

value within baseline period closest to and including day 0

### Follow up HbA1c:


value within follow up period closest to +180 days

## Results

**Study Cohort:** 1781 adults in the United States

Baseline Characteristics	
Age (years)	55 ± 10
Male	930 (52%)
Race	
African American	286 (16%)
Asian	35 (2%)
Caucasian	1279 (72%)
Other/unknown	181 (10%)
Insulin use	
Bolus insulin (with or without basal insulin)	675 (38%)
Basal insulin only	563 (32%)
No insulin	543 (30%)
GLP-1 use	
Dulaglutide	771 (43%)
Liraglutide	391 (22%)
Semaglutide	401 (23%)
Other GLP-1	218 (12%)
Baseline HbA1c (%)	9.8 ± 1.5
Median time from 1 <sup>st</sup> GLP-1 to FreeStyle Libre (days)	422 Inter-quartile range: 91-786

**Change in HbA1c:** After FSL acquisition, average HbA1c decreased from  $9.8 \pm 1.5\%$  to  $8.3 \pm 1.6\%$  (paired change of  $-1.5 \pm 2.0\%$ ,  $P < 0.001$ ). Significant improvements in HbA1c ( $P < 0.001$ ) were also observed in subgroups based on (A) insulin therapy, (B) GLP-1 type, and (C) duration between starting GLP-1 and FSL.



## Conclusions

In a large, real-world study of adults with T2D, those on prior GLP-1 therapy experienced significant, additional HbA1c improvements after acquiring FSL, irrespective of GLP-1 duration, GLP-1 type, or insulin therapy type. These findings support use of FSL in adults with T2D treated with GLP-1.